

Product Datasheet

Coronavirus (COVID-19 Spike Protein S2; 686-1208) Antigen, Recombinant >95% LA623

| | |
|----------------------------|---|
| Product Name | Coronavirus (COVID-19 Spike Protein S2; 686-1208) Antigen, Recombinant >95% |
| Catalog Number | LA623 |
| Tested Applications | LF, ELISA, CLIA |
| Brand | EastCoast Bio |
| Form/Appearance | Affinity Purified Recombinant Antigen (hFc-tag) |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Purity | >95% |
| Buffer | 200 mM Phosphate Buffer, pH 7.4 |
| Molecular Weight | 84 kDa |
| Expression System | Mammalian Cell (HEK293) |