

Product Datasheet

Coronavirus (COVID-19 Spike Protein S1; 1-674) Antigen, Recombinant 85% LA631

Product Name	Coronavirus (COVID-19 Spike Protein S1; 1-674) Antigen, Recombinant 85%
Catalog Number	LA631
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (His-tag)
Storage	Short Term (<1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing
Note	Peptide composed of amino acids 1-674 of spike protein S1 domain.
Purity	>85%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Molecular Weight	100 kDa
Expression System	Insect Cell
Preservatives	0.05% Sodium Azide