

Product Datasheet

Coronavirus (COVID-19 Spike Protein S1; 1-674) Antigen, Recombinant 85% LA634

Product Name Coronavirus (COVID-19 Spike Protein S1; 1-674) Antigen,
Recombinant 85%

Catalog Number LA634

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Affinity Purified Recombinant Antigen (His-tag)

Storage Short Term (<1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing

Note Peptide composed of amino acids 1-674 of spike protein S1 domain.

Purity >85%

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Molecular Weight 100 kDa

Expression System Mammalian Cell (HEK293)

Preservatives 0.05% Sodium Azide