

Product Datasheet

Coronavirus (COVID-19 Spike Protein; Full Length) Antigen, Recombinant >90% LA636

Product Name	Coronavirus (COVID-19 Spike Protein; Full Length) Antigen, Recombinant >90%
Catalog Number	LA636
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Full Length Recombinant Antigen
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of full length spike protein sequence.
Purity	≥90%
Buffer	Tris Buffer with salts and additives, pH 8.0
Molecular Weight	138 kDa
Expression System	Mammalian Cell (HEK293)