

Product Datasheet

Coronavirus (CoV-NP 229E; Full Length) Antigen, Recombinant >90% LA638

Product Name	Coronavirus (CoV-NP 229E; Full Length) Antigen, Recombinant >90%
Catalog Number	LA638
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of full length Coronavirus 229E nucleoprotein sequence.
Purity	>90%
Buffer	Tris Buffer with salts and additives, pH 8.0
Molecular Weight	47 kDa
Expression System	E. coli