

Product Datasheet

Angiotensin-converting Enzyme 2 (ACE2) Antigen, Recombinant ≥90% LA649

Product Name	Angiotensin-converting Enzyme 2 (ACE2) Antigen, Recombinant ≥90%
Catalog Number	LA649
Tested Applications	LF, ELISA, CMIA
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (His-tag)
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Purity	≥90%
Buffer	10 mM Phosphate Buffer, 150 mM Sodium Chloride, pH 7.4
Molecular Weight	90 kDa
Expression System	Mammalian Cell (CHO)