

## Product Datasheet

### Angiotensin-converting Enzyme 2 (ACE2) Antigen, Recombinant ≥90% LA650

<b>Product Name</b>	Angiotensin-converting Enzyme 2 (ACE2) Antigen, Recombinant ≥90%
<b>Catalog Number</b>	LA650
<b>Tested Applications</b>	LF, ELISA, CMIA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Affinity Purified Recombinant Antigen (hFc-tag)
<b>Storage</b>	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	≥90%
<b>Buffer</b>	10 mM Phosphate Buffer, 150 mM Sodium Chloride, pH 7.4
<b>Molecular Weight</b>	250 kDa
<b>Expression System</b>	Mammalian Cell (CHO)