

Product Datasheet

Progesterone-3-CMO-BSA Antigen LA651

Product Name	Progesterone-3-CMO-BSA Antigen
Catalog Number	LA651
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservatives	0.05% Sodium Azide