

## Product Datasheet

### Estrone-3-Glucuronide (E1G)-BSA Antigen LA655

|                            |   |
|----------------------------|---|
| <b>Product Name</b>        | Estrone-3-Glucuronide (E1G)-BSA Antigen                                   |
| <b>Catalog Number</b>      | LA655   |
| <b>Tested Applications</b> | LF, ELISA   |
| <b>Brand</b>               | EastCoast Bio   |
| <b>Form/Appearance</b>     | Liquid; Antigen-BSA conjugate   |
| <b>Storage</b>             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 7.4                                   |
| <b>Matched Pairs</b>       | HM1148, HM1149 or HM1150 (detection)                                      |
| <b>Preservatives</b>       | 0.05% Sodium Azide  |