

## Product Datasheet

### Risperidone-BSA Antigen LA656

<b>Product Name</b>	Risperidone-BSA Antigen
<b>Catalog Number</b>	LA656
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Storage</b>	-20°C
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Matched Pairs</b>	HM1151, HM1152, HM1153 or HM1154 (detection)
<b>Preservatives</b>	0.01% Sodium Azide