

Product Datasheet

Coronavirus (CoV-NP NL63) Antigen, Recombinant >95% LA657

Product Name	Coronavirus (CoV-NP NL63) Antigen, Recombinant >95%
Catalog Number	LA657
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen
Storage	Short Term (\leq 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of 130 total amino acids from C-terminus.
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 with 25 mM Potassium Carbonate
Expression System	E. coli