

Product Datasheet

Coronavirus (CoV-NP HKU1; Full Length) Antigen, Recombinant >90% LA658

Product Name	Coronavirus (CoV-NP HKU1; Full Length) Antigen, Recombinant >90%
Catalog Number	LA658
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (His-tag)
Storage	Short Term (\leq 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of 130 total amino acids from C-terminus.
Purity	>90% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 with 25 mM Potassium Carbonate
Molecular Weight	50 kDa
Expression System	E. coli