

Product Datasheet

Coronavirus (CoV-NP OC43; Full Length) Antigen, Recombinant >95% LA659

Product Name	Coronavirus (CoV-NP OC43; Full Length) Antigen, Recombinant >95%
Catalog Number	LA659
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (His-tag)
Storage	Short Term (\leq 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of amino acids 31-417 (no signal sequence)
Purity	>90% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 with 25 mM Potassium Carbonate
Molecular Weight	50 kDa
Expression System	E. coli