

Product Datasheet

Coronavirus (MERS-CoV Spike Protein S1; 367-606) Antigen, Recombinant >95% LA660

Product Name	Coronavirus (MERS-CoV Spike Protein S1; 367-606) Antigen, Recombinant >95%
Catalog Number	LA660
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen (His-tag)
Storage	Short Term (\leq 1 week): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Peptide composed of amino acids 367-606 of spike protein S1 domain.
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 with 25 mM Arginine
Expression System	E. coli
Preservatives	0.05% Sodium Azide