

Product Datasheet

Synuclein, alpha (α -Syn) Antigen, Recombinant >85% LA662

Product Name	Synuclein, alpha (α -Syn) Antigen, Recombinant >85%
Catalog Number	LA662
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Recombinant Antigen
Storage	$\leq -20^{\circ}\text{C}$. Avoid repeated freezing and thawing.
Purity	>85% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2 with 20% glycerol
Matched Pairs	HM1089, HM1092, HM1093 or HM1094
Expression System	E. coli
Preservatives	0.05% Sodium Azide