

## Product Datasheet

### Neurofilament L (NfL; 1-544) Antigen, Recombinant >85% LA665

<b>Product Name</b>	Neurofilament L (NfL; 1-544) Antigen, Recombinant >85%
<b>Catalog Number</b>	LA665
<b>Tested Applications</b>	LF, ELISA
<b>Alternative Names</b>	Alzheimer's disease (AD), Traumatic brain injury (TBI), Parkinson's disease (PD), Huntington's disease, Multiple sclerosis (MS). Amyotrophic lateral sclerosis (ALS)
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Affinity Purified Recombinant Antigen (His-tag)
<b>Storage</b>	≤ -20°C. Avoid repeated freezing and thawing.
<b>Purity</b>	>85%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2 with 20% glycerol
<b>Molecular Weight</b>	63 kDa
<b>Matched Pairs</b>	HM1095, HM1096, HM1097, HM1099, HM1100, HM1101
<b>Expression System</b>	Insect Cell
<b>Preservatives</b>	0.05% Sodium Azide