

Product Datasheet

Neurofilament L (NfL; 1-544) Antigen, Recombinant >85% LA665

| | |
|----------------------------|---|
| Product Name | Neurofilament L (NfL; 1-544) Antigen, Recombinant >85% |
| Catalog Number | LA665 |
| Tested Applications | LF, ELISA |
| Alternative Names | NfL |
| Brand | EastCoast Bio |
| Form/Appearance | Affinity Purified Recombinant Antigen (His-tag) |
| Storage | ≤ -20°C. Avoid repeated freezing and thawing. |
| Purity | >85% |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.2 with 20% glycerol |
| Molecular Weight | 63 kDa |
| Matched Pairs | HM1095, HM1096, HM1097, HM1099, HM1100, HM1101 |
| Expression System | Insect Cell |
| Preservatives | 0.05% Sodium Azide |