

Product Datasheet

Estrone-3-Glucuronide (E1G)-BSA Antigen LA668

Product Name	Estrone-3-Glucuronide (E1G)-BSA Antigen
Catalog Number	LA668
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM1148, HM1149 or HM1150 (detection)
Preservatives	0.05% Sodium Azide