

## Product Datasheet

### Estrone-3-Glucuronide (E1G)-BSA Antigen LA669

<b>Product Name</b>	Estrone-3-Glucuronide (E1G)-BSA Antigen
<b>Catalog Number</b>	LA669
<b>Tested Applications</b>	LF, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Antigen-BSA conjugate
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Matched Pairs</b>	HM1148, HM1149 or HM1150 (detection)
<b>Preservatives</b>	0.05% Sodium Azide