

Product Datasheet

Coronavirus (COVID-19 UK (Alpha) Variant B.1.1.7 Spike Protein; Full Length) Antigen, Recombinant >95% LA770

Product Name	Coronavirus (COVID-19 UK (Alpha) Variant B.1.1.7 Spike Protein; Full Length) Antigen, Recombinant >95%
Catalog Number	LA770
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Full Length Recombinant Antigen (His-tag, Strep Tag)
Storage	<-70°C. Avoid repeated freezing and thawing
Note	Peptide composed of full length spike protein variant B.1.1.7 (UK variant) sequence.
Purity	>95%
Buffer	PBS
Molecular Weight	180 kDa
Expression System	Mammalian Cell (HEK293)
Preservatives	0.05% Sodium Azide