

Product Datasheet

Coronavirus (COVID-19 Brazil Variant (Gamma) P.1/B.1.1.28 Spike Protein; Full Length) Antigen, Recombinant >95% LA773

Product Name	Coronavirus (COVID-19 Brazil Variant (Gamma) P.1/B.1.1.28 Spike Protein; Full Length) Antigen, Recombinant >95%
Catalog Number	LA773
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Full Length Recombinant Antigen (His-tag, Strep Tag)
Storage	<-70°C. Avoid repeated freezing and thawing
Note	Peptide composed of full length spike protein variant P.1/B.1.28 (Brazil variant) sequence.
Purity	>95%
Buffer	PBS
Molecular Weight	170 kDa.
Expression System	Mammalian Cell (HEK293)