

Product Datasheet

Diphtheria Toxoid Antigen, Native >98% LA776

Product Name	Diphtheria Toxoid Antigen, Native >98%
Catalog Number	LA776
Tested Applications	Control
Brand	EastCoast Bio
Form/Appearance	Liquid; Native Antigen
Storage	Short Term (≤ 1 weeks): 2-8°C. Long Term: -70°C. Avoid repeated freezing and thawing.
Note	No cytotoxicity or colony growth observed.
Purity	>98%
Source	Corynebacterium diphtheria toxoid from culture
Buffer	10 mM Tris, pH 7.5 with EDTA