

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-gp41) Antigen, Recombinant >90% LA779

Product Name	Human Immunodeficiency Virus 1 (HIV 1-gp41) Antigen, Recombinant >90%
Catalog Number	LA779
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen (His-tag)
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Epitope	gp41
Purity	>90%
Buffer	20 mM Tris-HCl, pH 8.0 with 1 mM EDTA
Specificity	gp41
Molecular Weight	MW= 15.8 kDa
Matched Pairs	LA038 (detection)
Expression System	E. coli