

Product Datasheet

Human Immunodeficiency Virus 1 (HIV 1-gp41) Antigen, Recombinant >85% LA780

Product Name	Human Immunodeficiency Virus 1 (HIV 1-gp41) Antigen, Recombinant >85%
Catalog Number	LA780
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	2-8°C
Epitope	gp41
Purity	>85%
Buffer	20 mM Carbonate-bicarbonate Buffer, pH 9.5
Specificity	gp41
Molecular Weight	MW= 32 kDa
Matched Pairs	LA782 (detection)
Expression System	E. coli
Preservatives	0.1% Sodium Azide