

Product Datasheet

Human Immunodeficiency Virus 2 (HIV 2-gp36) Antigen, Recombinant >80% LA781

Product Name	Human Immunodeficiency Virus 2 (HIV 2-gp36) Antigen, Recombinant >80%
Catalog Number	LA781
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Recombinant Antigen
Storage	2-8°C
Epitope	gp36
Purity	>80%
Buffer	20 mM Carbonate-bicarbonate Buffer, pH 9.5
Specificity	gp36
Molecular Weight	MW= 35 kDa
Matched Pairs	LA782 (detection)
Expression System	E. coli
Preservatives	0.1% Sodium Azide