

## Product Datasheet

### Human Immunodeficiency Virus 1 & 2 (HIV 1 & 2-gp36, gp41, gp120) Antigen, Recombinant >90% LA782

<b>Product Name</b>	Human Immunodeficiency Virus 1 & 2 (HIV 1 & 2-gp36, gp41, gp120) Antigen, Recombinant >90%
<b>Catalog Number</b>	LA782
<b>Tested Applications</b>	LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Recombinant Antigen
<b>Storage</b>	2-8°C
<b>Epitope</b>	gp36, gp41, gp120
<b>Purity</b>	>90%
<b>Buffer</b>	20 mM Tris, pH 8.5 with 5 mM EDTA
<b>Specificity</b>	gp36, gp41, gp120
<b>Molecular Weight</b>	MW= 43 kDa
<b>Matched Pairs</b>	LA780 or LA781 (capture)
<b>Expression System</b>	E. coli