

Product Datasheet

Lysergic acid diethylamide (LSD)- BSA Antigen LA783

Product Name	Lysergic acid diethylamide (LSD)- BSA Antigen
Catalog Number	LA783
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Antigen-BSA conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Buffer	PBS, pH 7.2
Matched Pairs	HM1201 (detection)
Preservatives	0.02% Sodium Azide