

Product Datasheet

Treponema pallidum (TpN15) Antigen, Recombinant >90% LA785

Product Name Treponema pallidum (TpN15) Antigen, Recombinant >90%

Catalog Number LA785

Tested Applications ELISA, CMIA

Alternative Names Syphilis

Brand EastCoast Bio

Form/Appearance Purified Recombinant Antigen

Storage -20°C. Avoid repeated freezing and thawing.

Purity ≥95% by SDS-PAGE

Buffer 20 mM Tris-HCl, pH 8.0 with 100 mM NaCl and 1 mM EDTA

Specificity TpN15

Expression System E. coli