

## Product Datasheet

### Hepatitis B Surface Antigen (HBsAg-subtype Ad), Native >95% LA787

**Product Name** Hepatitis B Surface Antigen (HBsAg-subtype Ad), Native >95%

**Catalog Number** LA787

**Tested Applications** ELISA

**Brand** EastCoast Bio

**Form/Appearance** Purified & 0.2 µm filtered Native Antigen

**Storage** 2-8°C

**Note** Tested negative for HIV 1/2 and HCV, however as with all biological material it should be handled as potentially infectious.

**Purity** >95%

**Source** Human Plasma

**Buffer** 25 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl, 5% sucrose

**Preservatives** 0.1% Sodium Azide