

Product Datasheet

HighScript Reverse Transcriptase #6601

Product Name HighScript Reverse Transcriptase

Catalog Number #6601

Description HighScript Reverse Transcriptase is a modified version of MMLV reverse transcriptase with noticeable thermostability and high enzymatic activity. This enzyme is offered as a blend with an RNase inhibitor to prevent RNA degradation. HighScript Reverse Transcriptase together with enhanced buffer chemistry enables fast, efficient and accurate synthesis of the cDNA molecule.

Tested Applications End-Point, Real-Time

Storage -20°C

Content 10 000 U: 200 U/μl, 2 x 0.025 ml HighScript Reverse Transcriptase; 1 x 0.2 ml 5x HighScript buffer
40 000 U: 200 U/μl, 2 x 0.1 ml HighScript Reverse Transcriptase; 4 x 0.2 ml 5x HighScript buffer