

## Product Datasheet

### qProbe Mix Separate ROX #2401

<b>Product Name</b>	qProbe Mix Separate ROX
<b>Catalog Number</b>	#2401
<b>Description</b>	qProbe Mix is a universal one-step probe mix for robust, sensitive, and fast qPCR. The mix uses state-of-the-art technologies with an antibody-regulated hot-start Taq DNA polymerase for real-time PCR amplification of single or multiplex DNA targets. qProbe Mix is compatible with several probes such as TaqMan® and Scorpions®.
<b>Tested Applications</b>	End-Point, Real-Time
<b>Brand</b>	Medix Biochemica
<b>Storage</b>	-20°C
<b>Content</b>	100 rxn: 1 x 1 ml qProbe Mix No-ROX (2x); 1 x 0.2 ml 50µM ROX Additive 500 rxn: 5 x 1 ml qProbe Mix No-ROX (2x); 1 x 0.2 ml 50µM ROX Additive