

## Product Datasheet

### qProbe Mix Separate ROX #2401

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | qProbe Mix Separate ROX  |
| <b>Catalog Number</b>      | #2401  |
| <b>Description</b>         | qProbe Mix is a universal one-step probe mix for robust, sensitive, and fast qPCR. The mix uses state-of-the-art technologies with an antibody-regulated hot-start Taq DNA polymerase for real-time PCR amplification of single or multiplex DNA targets. qProbe Mix is compatible with several probes such as TaqMan® and Scorpions®. |
| <b>Tested Applications</b> | End-Point, Real-Time   |
| <b>Brand</b>               | Medix Biochemica   |
| <b>Storage</b>             | -20°C  |
| <b>Content</b>             | 100 rxn: 1 x 1 ml qProbe Mix No-ROX (2x); 1 x 0.2 ml 50µM ROX Additive<br>500 rxn: 5 x 1 ml qProbe Mix No-ROX (2x); 1 x 0.2 ml 50µM ROX Additive   |