

Product Datasheet

Anti-Human CD25 Monoclonal Antibody, Azide Free Clone B-F2 852.010

Product Name	Anti-Human CD25 Monoclonal Antibody, Azide Free Clone B-F2
Catalog Number	852.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, Functional assay, IHC, FC
Immunogen	PHA activated T cells
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	PHA activated T cells
Specificity	Recognises the Interleukin 2 Receptor α chain (IL-2 Ra), a 55 kDa protein
Biological Activity	Inhibits mixed lymphocyte culture
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells