

Product Datasheet

qPCR Lyo-Ready Mix #4401

Product Name qPCR Lyo-Ready Mix

Catalog Number #4401

Description qPCR Lyo-Ready Mix is a universal probe mix that allows robust, sensitive, and fast qPCR for the detection of single or multiplex DNA targets. qPCR Lyo-Ready Mix contains optimized excipients, buffer, PCR enhancers and an antibody-regulated hot-start Taq DNA polymerase and is ready to be lyophilized to produce stable reagents at room temperature. Upon addition of target specific primers/probes to the mastermix, the mixture can be lyophilized directly, without the need to add additional excipients. qPCR Lyo-Ready Mix is compatible with several probes such as TaqMan® and Scorpions®.

Tested Applications End-Point, Real-Time, Lyophilized

Brand Medix Biochemica

Storage -20°C

Content 500 rxn: 4 x 1.25 ml qPCR Lyo-Ready Mix (2x)
10 000 rxn: 100 ml qPCR Lyo-Ready Mix (2x)