

Product Datasheet

Luteinizing Hormone (LH) Antibody OO44

Product Name	Luteinizing Hormone (LH) Antibody
Catalog Number	OO44
Species/Host	Mouse
Tested Applications	ELISA
Immunogen	Human pituitary LH
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate, pH 7.4 with 150 mM NaCl
Specificity	LH
Preservatives	0.05% Sodium Azide