

Product Datasheet

Anti-Human CD25 Monoclonal Antibody, Unconjugated Clone B-F2 852.013.020

Product Name Anti-Human CD25 Monoclonal Antibody, Unconjugated Clone B-F2

Catalog Number 852.013.020

Description Mouse monoclonal

Species/Host Mouse

Conjugation Purified

Tested Applications IHC, FC

Immunogen PHA activated T cells

Brand Diaclone

Form/Appearance Phosphate-buffered saline with 1% BSA and 0.09% sodium azide

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Note Use 10 µl to label 10⁶ cells or 100 µl of whole blood

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Source PHA activated T cells

Specificity Recognises the Interleukin 2 Receptor α chain (IL-2 Ra), a 55 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells