

Product Datasheet

Methamphetamine Antibody P01-99-31R-IF

Product Name	Methamphetamine Antibody
Catalog Number	P01-99-31R-IF
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Methamphetamine-BgG
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Methamphetamine-BgG linked at the Para Position off the Phenyl Group.
Clonality	Polyclonal
Buffer	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Cross Reactivity	Low cross reactivity to Amphetamine, Phenylpropanolamine
Specificity	Methamphetamine
Matched Pairs	P81-99-31A (capture) or P91-99-31H (detection)
Preservatives	0.05% Sodium Azide