

## Product Datasheet

### Fentanyl Antibody P01-99-53R-IF

<b>Product Name</b>	Fentanyl Antibody
<b>Catalog Number</b>	P01-99-53R-IF
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Fentanyl-BgG
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Polyclonal Antibody, IgG Fraction
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Fentanyl linked off the aliphatic group, ethyl. Additional lots available, please inquire.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	0.05 M Phosphate, pH 7.15 with 0.85% Sodium Chloride
<b>Specificity</b>	Fentanyl
<b>Matched Pairs</b>	P81-99-51A (capture) or P91-99-51H (detection)
<b>Preservatives</b>	0.05% Sodium Azide