

Product Datasheet

Methaqualone Antibody P02-99-83R-IF

Product Name	Methaqualone Antibody
Catalog Number	P02-99-83R-IF
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Methaqualone-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Specificity	Methaqualone
Matched Pairs	P82-99-81A (capture)
Preservatives	0.05% Sodium Azide