

Product Datasheet

Diphenhydramine Antibody P02-99-92R-IF

Product Name	Diphenhydramine Antibody
Catalog Number	P02-99-92R-IF
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Diphenhydramine -BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05 M Phosphate 0.85% Sodium Chloride, pH 7.15
Specificity	Diphenhydramine
Matched Pairs	P82-99-91A (capture)
Preservatives	0.05% Sodium Azide