

Product Datasheet

Herpes Simplex Virus (HSV-ICP4) Antibody P1101

Product Name	Herpes Simplex Virus (HSV-ICP4) Antibody
Catalog Number	P1101
Species/Host	Mouse
Tested Applications	WB, IFA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	≤-20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	ICP4