

Product Datasheet

Herpes Simplex Virus (HSV-gD) Antibody P1103

Product Name	Herpes Simplex Virus (HSV-gD) Antibody
Catalog Number	P1103
Species/Host	Mouse
Tested Applications	WB, IFA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	≤-20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	gD