

Product Datasheet

Anti-Human CD126 Monoclonal Antibody, Azide Free Clone B-F19 852.040

Product Name Anti-Human CD126 Monoclonal Antibody, Azide Free Clone B-F19

Catalog Number 852.040, 852.040.000, 852.040.005

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, IHC, FC

Immunogen gp80 transfected CHO cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm.
Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the Interleukin-6 Ra subunit, IL-6 Ra, IL-6 R, a 80 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells