

Product Datasheet

Herpes Simplex Virus 2 (HSV 2-gC) Antibody P1124

Product Name Herpes Simplex Virus 2 (HSV 2-gC) Antibody

Catalog Number P1124

Species/Host Mouse

Tested Applications IFA

Brand EastCoast Bio

Form/Appearance Protein G Purified Monoclonal Antibody

Storage $\leq -20^{\circ}\text{C}$. Avoid repeated freezing and thawing.

Isotype IgG2a, kappa

Clonality Monoclonal

Buffer 10 mM Phosphate Buffered Saline, pH 7.4

Specificity gC-2