

Product Datasheet

Cytomegalovirus (CMV-IE1) Antibody P1203

Product Name	Cytomegalovirus (CMV-IE1) Antibody
Catalog Number	P1203
Species/Host	Mouse
Tested Applications	WB, IFA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	≤-20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	IE1