

Product Datasheet

Cytomegalovirus (CMV-UL83) Antibody P1205

Product Name	Cytomegalovirus (CMV-UL83) Antibody
Catalog Number	P1205
Species/Host	Mouse
Tested Applications	WB, IFA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	≤-20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	UL83