

## Product Datasheet

### Cytomegalovirus (CMV-pp28) Antibody P1207

<b>Product Name</b>	Cytomegalovirus (CMV-pp28) Antibody
<b>Catalog Number</b>	P1207
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	WB, IFA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein G Purified Monoclonal Antibody
<b>Storage</b>	≤-20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a, kappa
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Specificity</b>	pp28