

Product Datasheet

Cytomegalovirus (CMV-IE1/2) Antibody P1215

Product Name	Cytomegalovirus (CMV-IE1/2) Antibody
Catalog Number	P1215
Species/Host	Mouse
Tested Applications	WB, IFA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	≤-20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	IE1/2