

Product Datasheet

Advanced Glycosylation End Product (AGE) Antibody PP102IF

Product Name	Advanced Glycosylation End Product (AGE) Antibody
Catalog Number	PP102IF
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Bovine Serum Albumin/Human Serum Albumin - Advanced Glycosylation Endproduct
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05M Phosphate, 0.85% NaCl, @ pH 7.15
Cross Reactivity	BSA-AGE; HAS-AGE (less than 1%)
Preservatives	0.05% Sodium Azide