

Product Datasheet

Cannabinoid (THC) Antibody PP124IF

Product Name	Cannabinoid (THC) Antibody
Catalog Number	PP124IF
Species/Host	Sheep
Tested Applications	LF, ELISA
Immunogen	Δ 8 THC-BgG
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody
Storage	Short Term (\leq 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	50 mM Phosphate, pH 7.15
Matched Pairs	LA020 (capture)
Preservatives	0.02% Sodium Azide